

BURNAID

Sterile Gel Dressing


- Ideal for larger burns
 - Suitable for transport
 - Sterile
 - Provides a physical barrier against contamination
 - Can be placed over burnt clothes
-
- 4 sizes – 10cmx10cm, 20cmx20cm, 55cmx40cm, Face Mask 30x40cm



Locate the nick at the top right hand front corner of the sachet. Place pinched thumb and forefingers either side of the 'nick' and pull both hands apart at 180 degrees to each other. Pull a strip down along the top of the dressing. Alternatively cut open with scissors



10x10 cm



20x20 cm



55x40 cm



40x30 cm face mask



Fast effective cooling and pain relief? The choice is natural